HemoCue® Plasma/Low Hb

Operating manual

Manual de funcionamiento

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 **HEMOCUE®**

GB/US

Table of Contents

HemoCue Plasma/Low Hb system	5
Components.....	6
Start up	8
Measuring Sample material	10
Maintenance	14
Troubleshooting Guide.....	16
Specifications.....	25

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Système HemoCue® Plasma/Low Hb

GB/US

Thank you for choosing the HemoCue® Plasma/Low Hb system. The system is used for quantitative determination of low levels of hemoglobin in plasma and serum specimens, aqueous solutions, or stored or banked erythrocytes using a specially designed photometer, the HemoCue Plasma/Low Hb Photometer and specially designed microcuvettes, the HemoCue Plasma/Low Hb Microcuvettes.

The HemoCue Plasma/Low Hb Photometer is only to be used with HemoCue Plasma/Low Hb Microcuvettes. The HemoCue Plasma/Low Hb Microcuvettes are for *In-Vitro* diagnostic use only.

All system components are designed and manufactured to provide maximum safety. Any other use of the system may impair the safety.

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Components
Componentes
Composants

GB/US

1. HemoCue Plasma/Low Hb Photometer*
2. AC adapter**
3. 5 alkaline type AA batteries***
4. A vial of HemoCue Plasma/Low Hb Microcuvettes***
5. HemoCue Plasma/Low Hb Operating Manual
6. HemoCue Cleaner

Open the carton and lift out the photo-meter and accessories.

If no AC power is available, use 5 alkaline type AA batteries. On the bottom of the photometer there is a lid covering the battery compartment. Press the flap to remove the lid. Place the batteries in the battery compartment and replace the lid.

Consult local environmental authorities for proper disposal of batteries.

* Do not open the cover of the photometer. Note: The warranty is void if the cover of the photometer has been opened.

**  Only use AC adapters listed under specifications.

*** Not included.

For information about HemoCue Plasma/Low Hb Microcuvettes please contact your HemoCue distributor.

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



Start up
Activación
Démarrage

GB/US

Make sure the photometer is placed horizontally and on a stable surface.

Only use AC adapters approved by HemoCue, as listed under 'Technical Specifications'.

1. If AC power is available, use the AC adapter.
2. If no power is available, insert 5 alkaline type AA batteries.
3. Turn the switch to the 'ON' position. All symbols appear on the display, and after approximately 10 seconds the letter 'LHb' will be displayed.
4. Pull out the cuvette holder to its loading position. This point, which should not be exceeded, is easily established by paying attention to a distinct stop.
5. After approximately 15 seconds the indication 'READY' appears on the display together with three flashing dashes. The photometer is ready for measurement.

Turn the photometer off by switching power switch to the 'OFF' position.

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
GB/US

Note: Always handle blood specimens with care, as they might be infectious. Always wear protective gloves when handling blood specimens.

1. To perform a test, the cuvette holder should be in its loading position. The display will show three flashing dashes and 'READY'. Remove a HemoCue Plasma/Low Hb Microcuvette from the vial.
2. Place a drop of well mixed sample onto a hydrophobic surface using a suitable transfer device.
3. Fill the microcuvette in one continuous process. Do NOT refill.
4. Wipe off excess sample from the outside of the microcuvette with a clean, lint-free wipe. Do not touch the filling end of the microcuvette. If a second sample is to be taken, fill a new microcuvette from a new drop of sample. This should not be done until the measurement of the first sample is completed.
5. Look for air bubbles in the filled microcuvette. If present, discard it and fill a new microcuvette from a new drop of sample. Small bubbles around the edge can be ignored.

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

GB/US

Note: Visually turbid samples should be filtered before analysis using a 0.2 µm filter.

- Place the microcuvette into the cuvette holder and start measurement as soon as possible but no later than 60 seconds after filling the microcuvette by gently pushing the cuvette holder to its measuring position. Note: Slamming the cuvette holder into place with undue force will cause splashing of the sample material onto the optical surfaces.
- During measurement 'MEASURING' and three fixed dashes will be shown.
- After approximately 15–60 seconds, the hemoglobin value is displayed. The result will remain on the display as long as the cuvette holder is in its measuring position. Do not remeasure the microcuvette.
- Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal.

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

Maintenance
Mantenimiento
Maintenance

GB/US

The cuvette holder should be cleaned after each day of use. Make sure that the photometer is turned off and that the display is blank.

1. Pull the cuvette holder completely out. Clean the cuvette holder with alcohol (20–70 %) or mild detergent. The cuvette holder can also be autoclaved. The cover may be cleaned with alcohol (20–70 %) or a mild detergent.
2. If the optical parts becomes dirty, an error code will be displayed. Before cleaning the optical parts, the cuvette holder should be completely removed according to section 1 above. To clean the optical parts inside the analyzer, push a HemoCue Cleaner swab into the opening of optic unit, but only as far as until the white part of the cleaner swab is inside the optic unit (Do not push the cleaner further than this for optimal cleaning). Move the cleaner back and fourth 5–10 times. If the cleaner swab becomes stained, repeat with a new cleaner. Wait 15 minutes before replacing the cuvette holder.

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Troubleshooting Guide

If you are unable to resolve the problem by following this Troubleshooting Guide, please contact your local HemoCue distributor or HemoCue AB. The photometer should be cleaned as recommended under 'Maintenance' prior to service or disposal. Consult local environmental authorities for proper disposal. The photometer has no serviceable parts. Do not open the cover of the photometer. Note: The warranty is void if the cover of the photometer has been opened.

Symptom	Explanation	Action
The photometer shows an error code.	May be a temporary fault.	Turn off the photometer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below.
ERROR 900	No stable endpoint of the measurement is found within the time range. 1. The microcuvette is faulty. 2. The circuit board is out of order.	1a. Check the expiration date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The photometer needs service. Contact your distributor.
ERROR 903	1. Disturbances on main power supply. 2. The optronic unit is out of order.	1 Change the wall socket or use battery power. 2. The photometer needs service. Contact your distributor.
ERROR 901–905	1. Dirty optical parts or faulty electronics or optronic unit.	1a. Turn off the photometer and clean the optical parts as described in the <i>Maintenance</i> section. 1b. The photometer needs service. Contact your distributor.
ERROR 906	1. Unstable blank value. The photometer might be cold.	1. Turn off the photometer and allow it to reach room temperature. If the problem continues, the photometer needs service. Contact your distributor.
ERROR 907	1. The battery power is too low.	1a. The batteries need to be replaced. Turn off the photometer and replace the batteries, 5 alkaline type AA. 1b. Use the AC adapter.
ERROR 908	The absorbance is too high. 1. Light blocking item in the cuvette holder	1a. Check that the photometer and microcuvettes are used according to the HemoCue Plasma/Low Hb operating manual and package insert. 1b. The photometer needs service. Contact your distributor.
ERROR 916	High turbidity sample	Filter visually turbid sample before analysis.
ERROR 918	Hardware or internal error	The photometer needs service. Contact your distributor.

Symptom	Explanation	Action
HHH	1. Measured value exceeds 30.0 g/L (3.00 g/dL, 3,000 mg/dL, 1.90 mmol/L).	1a. Check expiry dates of the microcuvettes. 1b. For sample above 30.0 g/L (3.00 g/dL, 3,000 mg/dL, 1.90 mmol/L) use the HemoCue 201+ system.
No characters on the display.	1. The photometer is not receiving power. 2. If on battery power, the batteries need to be replaced. 3. The display is out of order.	1a. Check that the AC adapter is connected to the photometer and the AC power supply. 1b. Check that the cable is not damaged. 2. Turn off the photometer and replace the batteries, 5 alkaline type AA. 3. The photometer needs service. Contact your distributor.
The display contains erroneous characters.	1. The display is out of order. 2. The microprocessor is out of order.	1. The photometer needs service. Contact your distributor. 2. The photometer needs service. Contact your distributor.
The display shows 'LOWBAT'.	1. The batteries need to be replaced. 2. If on AC power, the AC adapter or the circuit board is out of order.	1. Turn off the photometer and replace the batteries, 5 alkaline type AA. 2a. Check that the AC adapter is properly connected and working. 2b. The photometer needs service. Contact your distributor.
The display does not switch from 'LHb' to 'READY' or from 'READY' to 'MEASURING'.	1. The cuvette holder sensor is out of order.	1. The photometer needs service. Contact your distributor.

Symptom	Explanation	Action
Measurement on control materials are out of range — either too high or too low.	<ol style="list-style-type: none"> 1. The microcuvettes are beyond their expiry date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. The controls are beyond their expiry dates or have been improperly stored. 4. The control has not been mixed properly and/or is not at room temperature. 5. The microcuvette has not been placed in the photometer within 60 seconds of filling. 6. Air bubbles in the microcuvette. 7. The optical parts are dirty. 8. The control is not suitable for use with the HemoCue Plasma/Low Hb system. 9. The calibration of the photometer has been changed. 	<ol style="list-style-type: none"> 1. Check the expiry date and the storage conditions of the microcuvettes. 2. Remeasure the control with a new microcuvette. 3. Check the expiry date and the storage conditions of the control. Remeasure the control with a new microcuvette. If the problem continues, contact the manufacturer of the control. 4. Make sure that the control is mixed properly and at room temperature. If the problem continues, contact the manufacturer of the control. 5. Remeasure the control with a new microcuvette. 6. Check the microcuvette for air bubbles. Remeasure the control with a new microcuvette. 7. Clean the optical parts as described in the <i>Maintenance</i> section. 8. Contact your distributor for control information. 9. The photometer needs service. Contact your distributor.
Measurement on samples are higher or lower than anticipated.	<ol style="list-style-type: none"> 1. Improper sampling technique. 2. The microcuvettes are beyond their expiry date, damaged or have been improperly stored. 3. The optical eye of the microcuvette is contaminated. 4. The sample has not been mixed properly. 5. Air bubbles in the microcuvette. 6. The optical parts are dirty. 7. The calibration of the photometer has been changed. 	<ol style="list-style-type: none"> 1. See the <i>Measuring</i> section in this manual. 2. Check the expiry date and the storage conditions of the microcuvettes. 3. Remeasure the sample with a new microcuvette. 4. Make sure the sample is properly mixed. 5. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 6. Clean the optical parts as described in the <i>Maintenance</i> section. 7. The photometer needs service. Contact your distributor.

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Specifications

Intended Purpose/Intended Use

The HemoCue Plasma/Low Hb system is used for quantitative determination of low levels of hemoglobin in plasma and serum specimens, aqueous solutions, or stored or banked erythrocytes using a specially designed photometer, the HemoCue Plasma/Low Hb Photometer and specially designed microcuvettes, the HemoCue Plasma/Low Hb Microcuvettes. The HemoCue Plasma/Low Hb Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Plasma/Low Hb Photometer is only to be used with HemoCue Plasma/Low Hb Microcuvettes.

IVD Medical Device Directive

The HemoCue Plasma/Low Hb system complies with the IVD Medical Device Directive 98/79/EC and carry the CE mark.

Principle of the method/procedure

Principle of the method

The reaction in the microcuvette is a modified azidemethemoglobin reaction. The erythrocytes are hemolyzed to release the hemoglobin. The hemoglobin is converted to methemoglobin and then combined with azide to form azidemethemoglobin. The measurement takes place in the photometer in which the transmittance is measured and the absorbance and hemoglobin level is calculated. The absorbance is directly proportional to the hemoglobin concentration.

Principle of the procedure

The system consists of an photometer together with microcuvettes. The microcuvette serves both as a pipette and as a measuring cuvette and is for single-use only. A sample of approximately 20 µL is drawn into the cavity by capillary action. The photometer measures at two wavelengths in order to compensate for a certain degree of turbidity, and the hemoglobin level is calculated and presented (limitations apply as described in '*Limitation of the procedure*').

The HemoCue Plasma/Low Hb system is calibrated against the international reference method for hemoglobin determination, ICSH and needs no further calibration.

Warning and precautions

The microcuvettes are for *In Vitro* Diagnostic use only. Non-laboratory personnel should receive adequate training before using this system for the first time. Always handle blood specimens with care as they may be infectious. Consult local environmental authorities for proper disposal.

Storage and handling

Operating temperature for the system is 15–30°C (59–86°F). The system should not be operated at < 5 % or > 90 % non-condensing humidity.

HemoCue® Plasma/Low Hb Microcuvettes

The microcuvettes are to be stored at 15–30°C (59–86°F). Do not refrigerate. Use the microcuvettes prior to the expiration date that is printed on the package. Once the seal of the vial is broken, the microcuvettes are stable for three months. Keep the vial properly closed. All unused microcuvettes should remain in the original package.

HemoCue® Plasma/Low Hb Photometer

The photometer can be stored and transported in temperatures between 0–50°C (32–122°F). Allow the photometer to reach ambient temperature before use. The photometer should not be stored at < 5 % or > 90 % non-condensing humidity.

Specimen collection and preparation

Plasma/serum specimens and aqueous solutions containing hemoglobin, such as irrigating fluid from surgical procedures, may be used. Mix the specimen thoroughly before use. Supernatant from erythrocyte suspensions may be used. Carefully separate the supernatant from the erythrocytes using accepted laboratory procedure.

Materials required

- HemoCue Plasma/Low Hb Photometer
- HemoCue Plasma/Low Hb Microcuvettes
- Pipette or other transfer device
- Disposable pipette tips
- Lint-free wipe (non-fraying)
- Hydrophobic surface

Quality Control

Two levels of liquid control is recommended to be run on the day of use, or follow local guidelines regarding quality control procedures. Only use controls recommended by HemoCue, see relevant package insert for more information.

Measuring range

The system is linear from between 0.3–30.0 g/L (0.03–3.00 g/dL, 30–3,000 mg/dL, 0.02–1.90 mmol/L). Caution should be taken when evaluating instrument readings between 0 and 0.3 g/L (0–0.03 g/dL, 0–30 mg/dL, 0–0.02 mmol/L).

Results

The measured hemoglobin value is read directly from the HemoCue Plasma/Low Hb Photometer. No calculations are necessary. Applications sheets for use in determining blood loss and for calculating the amount of free hemoglobin in salvaged or stored blood are available from HemoCue AB.

Limitations of the procedure

- a) If 'HHH' is displayed, the result exceeds the measuring range of the system.
- b) For values above 30.0 g/L (3.00 g/dL, 3,000 mg/dL, 1.90 mmol/L) use a suitable laboratory method or the HemoCue Hb 201 system.
- c) Sulfhemoglobin is not measured with this method.
- d) Levels of bilirubin up to 340 µmol/L (20 mg/dL) do not influence the assay.
- e) The presence of lipid may interfere with the hemoglobin determination. Therefore, samples which are visibly turbid should be filtered (pore size 0.2 µm). If the user feels uncertain about the degree of turbidity, we recommend filtering the sample.
- f) The performance characteristics of this system have not been determined using samples obtained from uremic patients.

Specific performance characteristics

The results given below are from a HemoCue Plasma/Low Hb Photometer standardized against the HiCN-(ICSH)² method.

Within-run precision

Within-run precision was determined on 5 photometers and by using within house prepared aqueous and plasma solutions.

Hemoglobin in 0.9 % NaCl.

Solution	Mean (g/L)	SD (g/L)	CV %	Number
Low	1.59	0.06	3.6	20
Mean	16.30	0.14	0.8	20
High	25.30	0.16	0.6	20

Hemoglobin in heparinized plasma.

Solution	Mean (g/L)	SD (g/L)	CV %	Number
Low	1.95	0.08	4.2	20
Mean	15.30	0.11	0.7	20
High	25.50	0.22	0.9	20

Between-run precision

Between-run precision was determined by using within house prepared solutions. The sample were analyzed on 20 consecutive working days, twice a day, in duplicate.

Hemoglobin in 0.9 % NaCl.

Solution	Mean (g/L)	SD (g/L)	CV %	Number
Low	1.43	0.06	4.3	20
Mean	14.80	0.17	1.2	20
High	24.80	0.25	1.0	20

Hemoglobin in heparinized plasma.

Solution	Mean (g/L)	SD (g/L)	CV %	Number
Low	1.78	0.09	5.0	20
Mean	14.80	0.21	1.4	20
High	24.50	0.25	1.0	20

Correlation study

- 1) Aqueous samples from surgical procedures, analyzed on the HemoCue Plasma/Low Hb system with the HiCN-(ICSH) method as reference.
- 2) Erythrocyte suspensions, analyzed on the HemoCue Plasma/Low Hb system with the HiCN-(ICSH) method as reference.
- 3) Erythrocyte suspensions, analyzed on the HemoCue Plasma/Low Hb system with the HiCN-(ICSH) method as reference.

		Reference				HemoCue			
Group	N	Mean (g/L)	Min (g/L)	Max (g/L)	SD (g/L)	Mean (g/L)	Min (g/L)	Max (g/L)	SD (g/L)
1	42	4.04	0.94	15.18	0.047	4.05	0.84	15.64	0.059
2	20	1.47	0.65	2.49	0.066	1.66	0.70	2.75	0.130
3	64	2.54	0.40	6.45	0.122	2.65	0.40	6.70	0.084

The regression line and the coefficient of correlation

Group	N	Regression line	The coefficient of correlation	Range	Photometers
1	42	$1.039 \times \text{ICSH} - 0.150$	0.9997	0.94–15.18	5
2	20	$1.149 \times \text{ICSH} - 0.024$	0.9913	0.65–2.49	1
3	64	$1.015 \times \text{ICSH} - 0.078$	0.9981	0.40–6.45	1

Technical Specifications

Dimensions: 160 x 210 x 90 mm (6.29 x 8.26 x 3.54 inches)
Weight: 1,000 g (2.20 pounds)
5 alkaline type AA batteries
Power adapter: CE marked
Only use adapters, as listed under 'Adapters'.
Pollution degree: 2
Overvoltage category: II
Atmospheric pressure: 800 hPa to 1060 hPa.
Equipment not suitable for use in the presence of flammable mixtures.
The instrument is tested according to IEC 61010-1, Second edition: 2001 and EN 61010-1: 2001, IEC/EN 61010-2-101: 2002, IEC 61326-1: 2005 and EN 61326-1: 2006, IEC 61326-2-6, First edition: 2005 and EN 61326-2-6: 2006, IEC 60601-1-2, Third edition: 2007 and EN 60601-1-2: 2007 and complies with the IVD Medical Device Directive 98/79/EC.
The instrument is made for continuous mode.

Essential Performance

The essential performance is quantitative determination of hemoglobin in plasma and serum specimens, aqueous solutions, or stored banked erythrocytes.

AC adapters

Country: EU/US/GB
Type: 12 V Friwo FW7333SM/12 V HCA01
Input: 100 V~240 V~/50-60 Hz/ ≤ 200 mA

Recommended separation distance between Portable and mobile RF communications equipment and HemoCue® Plasma/Low Hb Photometer

The HemoCue systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of HemoCue systems can help prevent electromagnetic interference by maintaining a minimum distance between portable and RF communications equipment (transmitters)

and HemoCue systems as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at maximum output power not listed above, the recommended separation distances (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – Electromagnetic immunity

The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or user of the HemoCue systems should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	< 5 % U (> 95 % dip in U) for a 0.5 cycle 40 % U (60 % dip in U) for 5 cycles 70 % U (30 % dip in U) for 25 cycles < 5 % U (> 95 % dip in U) for 5 seconds For explanation of U see NOTE 1	< 5 % U (> 95 % dip in U) for a 0.5 cycle 40 % U (60 % dip in U) for 5 cycles 70 % U (30 % dip in U) for 25 cycles < 5 % U (> 95 % dip in U) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the HemoCue systems requires continued operation during power mains interruptions, it is recommended that the HemoCue systems be powered from an uninterruptible power supply or a battery.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
<p data-bbox="293 259 407 301">Conducted RF IEC 61000-4-6</p> <p data-bbox="293 326 407 368">Radiated RF IEC 61000-4-3</p>	<p data-bbox="532 259 678 301">3 Vrms 150 kHz to 80 MHz</p> <p data-bbox="532 326 678 368">3 V/m 80 MHz to 2.5 GHz</p> <p data-bbox="532 393 711 415">See NOTE 2 and NOTE 3</p>	<p data-bbox="769 259 824 281">3 Vrms</p> <p data-bbox="769 326 824 348">3 V/m</p>	<p data-bbox="943 259 1325 391">Portable and mobile RF communications equipment should be used no closer to any part of the HemoCue systems, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p data-bbox="943 393 1203 415">Recommended separation distance</p> <p data-bbox="943 441 1019 462">$d = 1.2\sqrt{P}$</p> <p data-bbox="943 488 1170 509">$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p data-bbox="943 535 1170 556">$d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p data-bbox="943 582 1325 665">Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in meters (m).</p> <p data-bbox="943 691 1325 774">Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).</p> <p data-bbox="943 800 1325 842">Interference may occur in the vicinity of equipment marked with the following symbol.</p> <div data-bbox="1101 867 1170 926" style="text-align: center;">  </div>

NOTE 1 U is the a.c. mains voltage prior to application of the test level.

NOTE 2 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 3 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HemoCue systems are used exceeds the applicable RF compliance level above, the HemoCue systems should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the systems.
- b) Over the frequency range 150 kHz to 80 MHz, field strength should be less than 3 V/m.

Technical specifications (EMC-RF)**Use only cables with the following specification:**

USB shielded maximum 2 m

Serial shielded maximum 1.5 m

Guidance and manufacturer's declaration — electromagnetic emissions		
The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or the user of the HemoCue systems should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment — guidance
RF emissions	Group 1	The HemoCue systems uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	The HemoCue systems are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

NOTE 1 It is the manufacturer's responsibility to provide equipment electromagnetic compatibility information to the customer or user.

NOTE 2 It is the user's responsibility to ensure that a compatible electromagnetic environment for the equipment can be maintained in order that the device will perform as intended.

Warranty

The photometer carries a 24-month warranty from the day of receipt. After the expiration date of the warranty, maintenance and repairs are offered at a fixed price. Any other use of the system than recommended by the manufacturer will void the warranty.

Service and disposal

The photometer should be cleaned as recommended under 'Maintenance' prior to service or disposal. Consult local environmental authorities for proper disposal.

Spare parts and accessories:

The following accessories and spare parts are available:

- Power adapter
- Battery lid
- Cuvette holder
- HemoCue Cleaners

Symbols used

Caution

CE mark

Class II equipment

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

Temperature limitation

Efficiency Level

Relative humidity, non-condensing

Consult instructions for use

To maintain safety use only adapter marked HCA01

DC inlet

'OFF' (power)

'ON' (power)

Biological risks

References

1. HemoCue Plasma/Low Hb package insert.
2. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A
3. ICSH Standard EP6/4: 1995, J Clin. Pathol. 1996; 49:271–274

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HemoCue AB | Box 1204 | SE-262 23 Ängelholm | Sweden | Phone +46 77 570 02 10 | Fax +46 77 570 02 12