IVD

150322 200320

GB/US

HemoCue® Plasma/Low Hb Microcuvettes

The HemoCue Plasma/Low Hb Microcuvettes are designed for use with the HemoCue Plasma/Low Hb Photometer. Please read the Operating Manual for proper use of the system¹.

Intended Purpose/Intended Use

The HemoCue Plasma/Low Hb system is used for quantitative determination of low levels of hemoglobin in plasma and serum specimens, aqueous solutions, or stored or banked erythrocytes using a specially designed photometer, the HemoCue Plasma/Low Hb Photometer and specially designed microcuvettes, the HemoCue Plasma/Low Hb Microcuvettes. The HemoCue Plasma/Low Hb Microcuvettes are for *In-Vitro* Diagnostic use only. The HemoCue Plasma/Low Hb Photometer is only to be used with HemoCue Plasma/Low Hb Microcuvettes.

IVD Medical Device Directive

The HemoCue Plasma/Low Hb Microcuvettes comply with the IVD Medical Device Directive 98/79/EC and carry the CE mark.

Principles of the method/procedure

Principle of the method

The reaction in the microcuvette is a modified azidemethemoglobin reaction. The erythrocytes are hemolyzed to release the hemoglobin. The hemoglobin is converted to methemoglobin and then combined with azide to form azidemethemoglobin. The measurement takes place in the photometer in which the transmittance is measured and the absorbance and hemoglobin level is calculated. The absorbance is directly proportional to the hemoglobin concentration.

Principle of the procedure

The system consists of a photometer together with microcuvettes. The microcuvette serves both as a pipette and as a measuring cuvette and is for single-use only. A sample of approximately 20 µL is drawn into the cavity by capillary action. The photometer measures at two wavelengths in order to compensate for a certain degree of turbidity, and the hemoglobin level is calculated and presented (limitations apply as described in "Limitation of the method"). The HemoCue Plasma/Low Hb system is calibrated against the international reference method for hemoglobin determination, ICSH² and needs no further calibration.

Composition

The microcuvette is made of polystyrene plastic. Reagents; <3000 µg/g microcuvette sodium deoxycholate, <1500 µg/g microcuvette sodium azide, <1650 µg/g microcuvette sodium nitrite, <700 µg/g microcuvette nonreactive ingredients.

Warning and precautions

The microcuvettes are for *In-Vitro* Diagnostic use only. Non-laboratory personnel should receive adequate training before using this system for the first time. Always handle blood specimens with care as they may be infectious. Consult local environmental authorities for proper disposal.

Storage and handling of the HemoCue Plasma/Low Hb Microcuvettes

The microcuvettes are to be stored at 15-30 °C (59-86 °F). Do not refrigerate. Use the microcuvettes prior to the expiration date that is printed on the package. Once the seal of the vial is broken, the microcuvettes are stable for three months. Keep the vial properly closed. All unused microcuvettes should remain in the original package.

Specimen collection and preparation

Plasma/serum specimens and aqueous solutions containing hemoglobin, such as irrigating fluid from surgical procedures, may be used. Mix the specimen thoroughly

before use. Supernatant from erythrocyte suspensions may be used. Carefully separate the supernatant from the erythrocytes using accepted laboratory procedure.


Directions for use

Materials required


- HemoCue Plasma/Low Hb Photometer
- HemoCue Plasma/Low Hb Microcuvettes
- Pipette or other transfer device
- Disposable pipette tips
- Lint-free wipe (non-fraying)
- Hydrophobic surface

Procedure


Please read the Operating Manual for proper use of the system¹.

1. Optical eye
2. Filling end

Dispense an aliquot of the well-mixed specimen onto a hydrophobic surface.

Remove the microcuvette from the package. Close the vial. Hold the microcuvette opposite the filling end and bring into contact with the specimen, allowing the cavity to fill completely. Always avoid touching the optical eye. Do not refill the cavity of the microcuvette. Wipe off the outside of the microcuvette with a clean lint-free wipe, being careful not to touch the open end. If air bubbles are seen in the optical eye of the microcuvette, discard the microcuvette and take a new sample. Small air bubbles around the edge do not influence the result. Place the microcuvette into the cuvette holder and start measurement as soon as possible but no later than 60 seconds after filling the microcuvette by gently pushing the cuvette holder to its measuring position. Within 60 seconds the result will appear on the display. Pull the cuvette holder out to its loading position and discard the used microcuvette.

Quality control

Two levels of liquid control is recommended to be run on the day of use, or as required by accrediting agencies or local policy. Contact HemoCue for current recommendations regarding controls.

Limitations of the method

- If "HHH" is displayed, the result exceeds the measuring range of the system.
- For values above 30.0 g/L (3.00 g/dL, 3000 mg/dL, 1.86 mmol/L) use a suitable laboratory method i.e. the HemoCue Hb 201 system.
- Sulfhemoglobin is not measured with this method.
- Levels of bilirubin up to 340 µmol/L (20 mg/dL) do not influence the assay.
- The presence of lipid may interfere with the hemoglobin determination. Therefore, samples which are visibly turbid should be filtered (pore size 0.2 µm). If the user feels uncertain about the degree of turbidity, we recommend filtering the sample.
- The performance characteristics of this system have not been determined using samples obtained from uremic patients.

Measuring range

The system is linear from between 0.3-30.0 g/L (0.03-3.00 g/dL, 30-3000 mg/dL, 0.02-1.86 mmol/L)³. Caution should be taken when evaluating instrument readings between 0 and 0.3 g/L (0-0.03 g/dL, 0-30 mg/dL, 0-0.02 mmol/L).

Results

The measured hemoglobin value is read directly from the HemoCue Plasma/Low Hb Photometer. No calculations are necessary. Applications sheets for use in determining blood loss and for calculating the amount of free hemoglobin in salvaged or stored blood are available from HemoCue AB.

Specific performance characteristics

The results given below are from a HemoCue Plasma/Low Hb Photometer standardized against the HiCN-(ICSH)² method.

Within-run precision

Within-run precision was determined on 5 photometers and by using within house prepared aqueous and plasma solutions.

Hemoglobin in 0.9 % NaCl.

Solution	Mean (g/L)	SD (g/L)	CV %	Number
Low	1.59	0.06	3.6	20
Mean	16.30	0.14	0.8	20
High	25.30	0.16	0.6	20

Hemoglobin in heparinized plasma.

Solution	Mean (g/L)	SD (g/L)	CV %	Number
Low	1.95	0.08	4.2	20
Mean	15.30	0.11	0.7	20
High	25.50	0.22	0.9	20

Between-run precision

Between-run precision was determined by using within house prepared solutions. The samples were analyzed on 20 consecutive working days, twice a day, in duplicate.

Hemoglobin in 0.9 % NaCl.

Solution	Mean (g/L)	SD (g/L)	CV %	Number
Low	1.43	0.06	4.3	20
Mean	14.80	0.17	1.2	20
High	24.80	0.25	1.0	20

Hemoglobin in heparinized plasma.

Solution	Mean (g/L)	SD (g/L)	CV %	Number
Low	1.78	0.09	5.0	20
Mean	14.80	0.21	1.4	20
High	24.50	0.25	1.0	20

Correlation study

- 1) Aqueous samples from surgical procedures, analyzed on the HemoCue Plasma/Low Hb system with the HiCN-(ICSH) method as reference.
- 2) Erythrocyte suspensions, analyzed on the HemoCue Plasma/Low Hb system with the HiCN-(ICSH) method as reference.
- 3) Erythrocyte suspensions, analyzed on the HemoCue Plasma/Low Hb system with the HiCN-(ICSH) method as reference.

		Reference				HemoCue			
Group	N	Mean (g/L)	Min (g/L)	Max (g/L)	SD (g/L)	Mean (g/L)	Min (g/L)	Max (g/L)	SD (g/L)
1	42	4.04	0.94	15.18	0.047	4.05	0.84	15.64	0.059
2	20	1.47	0.65	2.49	0.066	1.66	0.70	2.75	0.130
3	64	2.54	0.40	6.45	0.122	2.65	0.40	6.70	0.084

The regression line and the coefficient of correlation

Group	N	Regression Line	The coefficient of correlation	Range	Photometers
1	42	1.039xICSH-0.150	0.9997	0.94–15.18	5
2	20	1.149xICSH-0.024	0.9913	0.65– 2.49	1
3	64	1.015xICSH+0.078	0.9981	0.40– 6.45	1

Bibliography

1. HemoCue Plasma/Low Hb Operating Manual
2. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A
3. ICSH Standard EP6/4: 1995, J Clin. Pathol. 1996; 49:271-274

Symbols used

	Attention, see instructions for use		Batch code		Use by (Year Month Day)
	In vitro diagnostic medical device		Catalog number		Date of opening
	Do not reuse	15 °C 59 °F	30 °C 86 °F		Once opened – Expiry date
	CE mark	Temperature limitation			Open vial expiry date. Must not exceed the “Use by” date

PATENTED This product is covered by one or more patents

HemoCue Distributor USA

HemoCue America 250 South Kraemer Boulevard Brea, CA, 92821 USA	Phone (general) 800-881-1611 Orders 800-323-1674 Technical support 800-426-7256 Fax (cust. service) 800-333-7043 E-mail: info@hemocue.com www.hemocue.com
--	---

Manufacturer

HemoCue AB
Kuvettgatan 1
SE-262 71 Ängelholm
Sweden

Phone +46 77 570 0210
Fax +46 77 570 0212
E-mail info@hemocue.se
www.hemocue.com

HemoCue® Plasma/Low Hb Microcuvettes

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