150702 211012

US

HemoCue® Glucose 201 Microcuvettes

The HemoCue Glucose 201 Microcuvettes are designed for use with the HemoCue Glucose 201 Analyzer and the HemoCue Glucose 201 DM Analyzer (hereafter referred to as the HemoCue Glucose 201 Analyzer in this document). HemoCue Glucose 201 Microcuvettes are available in individual packages or in vials. Please read the relevant Operating Manual¹ for the analyzer for proper use of the system.

Intended use

Quantitative determination of glucose in whole blood using a specially designed analyzer, the HemoCue Glucose 201 Analyzer with plasma conversion. The quantitative determination of the instant blood glucose concentration in circulation supplements the clinical evidence in the diagnosis and treatment of patients with diabetes as well as monitoring of neonatal blood glucose levels. The HemoCue Glucose 201 Analyzer with plasma conversion multiplies the measured whole blood glucose value by a factor of 1.11² and displays a plasma equivalent glucose result. HemoCue Glucose 201 Microcuvettes are for *in vitro* diagnostic use only, and are only to be used with the HemoCue Glucose 201 Analyzer. **For professional use only.**

Summary and explanation of the test

As described in a review by Burrin and Alberti³, chemical micro methods for the determination of glucose in whole blood have been in use for decades. The oldest methods relied on chemistry based on the reducing properties of the 1,2-enediol form of glucose found under alkaline conditions. The next step in the development was to use chemistry methods based on ortho-toluidine. These chemistry methods were cumbersome and susceptible to interferences from a large number of chemical substances. Today's chemistry methods for the determination of glucose in body fluids are based on enzymes that are highly specific for glucose and ideally suited for use with photometric devices. Three different enzymes, glucose oxidase, hexokinase and glucose dehydrogenase can be used for the determination of glucose in body fluids. These enzymatic reactions can be adapted for use with a photometric device in the UV light range or in the visible light range after the addition of a chromogen to the reagents.

The disposable HemoCue Glucose 201 Microcuvette is made of plastic and contains a cavity that allows for drawing up approximately 5 μ L of specimen. The microcuvette contains dry reagents deposited uniformly within the cavity. The blood sample is drawn into the cavity by capillary action and is spontaneously mixed with the reagents. The chemistry method utilized by the HemoCue Glucose 201 Microcuvette is a modified glucose dehydrogenase method described by Banauch et al⁴. A chromogen compound is added to the reagents according to the principle outlined by Bergmeyer⁵, with saponin used for hemolyzing the erythrocytes. When used with the HemoCue Glucose 201 Analyzer, the end-point reaction within the microcuvette yields a colored compound proportional to the concentration of glucose present in the sample. The chemical end-point is reached, depending on the glucose concentration, in less than 90 seconds up to 155 mg/dL (8.6 mmol/L) glucose and less than 240 seconds at 444 mg/dL (24.6 mmol/L) glucose. The HemoCue Glucose 201 Microcuvettes are ready for use upon removal from the package. Capillary action draws 5 μ L of the sample into the cavity where the blood reacts with the reagents. The microcuvette is placed in the HemoCue Glucose 201 Analyzer and the chemical reaction is measured photometrically by the instrument until an end-point condition is reached. Utilization of a conversion factor allows the measured whole blood result to be displayed as a plasma equivalent glucose value. The HemoCue Glucose 201 Analyzer uses a dual wavelength measuring method at 667 nm and 840 nm. The HemoCue Glucose 201 Microcuvette is disposable once the test is completed.

Principles of the procedure

The technique

The HemoCue Glucose 201 technique is based on an optical measuring microcuvette of a small but precise volume and short light path. The microcuvette cavity contains reagents deposited on its inner walls. The blood sample is drawn into the cavity by capillary action and is mixed spontaneously with the reagents. The microcuvette is then placed in a HemoCue Glucose 201 Analyzer in which the transmittance is measured and the absorbance and plasma equivalent glucose level is calculated. Thus the technique makes it possible to sample the blood, mix and chemically react it with the reagents in the same microcuvette as is used for measurement.

The microcuvette

The microcuvette is made of polystyrene plastic and contains a cavity that holds approximately 5 μ L of specimen. The distance between the walls of the optical window is about 0.16 mm, which permits photometric determination of glucose in undiluted whole blood or control material.

The chemistry principle

The chemical reaction in the cavity has two phases; hemolysis, and the glucose reaction. Hemolysis, the disintegration of erythrocyte membranes, is caused by saponin. The glucose reaction is a modified glucose dehydrogenase method in which a tetrazolium salt is used to obtain a quantification of glucose in visible light. α -D-glucose is transformed to β -D-glucose using mutarotase. Glucose dehydrogenase acts as a catalyst for the oxidation of β -D-glucose, to form NADH, which in the presence of diaphorase produces a colored formazan with MTT, a tetrazolium salt. The formazan is quantified photometrically using a two wavelength photometric method at 667 nm and 840 nm.

Reagents

20 w/w % enzyme mix (mutarotase, glucose dehydrogenase, diaphorase, NAD), 20 w/w % MTT, 30 w/w % saponin, nonreactive ingredients.

Warning and precautions

HemoCue Glucose 201 Microcuvettes are for *in vitro* diagnostic use only. The chemicals deposited in the cavity of the microcuvette are harmful if swallowed. The microcuvettes should not be stored together with dry ice. Blood may be infectious, always handle with care and wear protective gloves. Follow local safety procedures for disposal of used microcuvettes.

Storage and handling of the HemoCue Glucose 201 Microcuvettes

The reagents within HemoCue Glucose 201 Microcuvettes are moisture and temperature sensitive. The color of the reagents in dry form is faint yellow. Improper storage may give microcuvettes a distinctive brown or blue-brown color and these microcuvettes should not be used. As the test method relies on photometric measurement, care should be taken not to hold the microcuvette by the filling end. Also take care to wipe away all contaminating substances from the outer surface of the microcuvette. All unused microcuvettes must remain in the original package. Use the microcuvettes prior to the expiry date. The expiry date is printed on each package.

Storage for microcuvettes kept in a vial

Store unopened HemoCue Glucose 201 Microcuvettes below 46 °F (8° C) (incl. storage in a freezer). Note, microcuvettes stored in a freezer must be allowed to reach room temperature 59-86 °F (15-30 °C) (approx. 30 minutes) before analysis. Microcuvettes kept in an opened vial are stable for 30 days when stored in a refrigerator at 35-46 °F (2-8 °C). After breaking the seal, the vial can be stored at room temperature 59-86 °F (15-30 °C) for up to 3 days. Close the lid immediately after microcuvettes are removed from the vial.

Storage for individually packaged microcuvettes

Store unopened HemoCue Glucose 201 Microcuvettes below 46 °F (8° C) (incl. storage in a freezer). Note, microcuvettes stored in a freezer must be allowed to reach room temperature 59-86 °F (15-30 °C) (approx. 30 minutes) before analysis. The individually packaged microcuvettes can be stored at room temperature 59-86 °F (15-30 °C) for up to 3 days.

Instrument

The HemoCue Glucose 201 Microcuvette is specifically designed for use with the HemoCue Glucose 201 Analyzer. The system is factory calibrated and needs no further calibration. The HemoCue glucose reference system is traceable to an Isotope Dilution Gas Chromatography - Mass Spectrometry (ID GC-MS) method.

Specimen collection and preparation

Capillary, venous or arterial whole blood may be used. The anticoagulants EDTA or sodium heparin and the glycolysis inhibitors sodium fluoride, sodium oxalate or potassium oxalate shall be used.

Glycolysis is a major concern in all glucose measurements. To minimize the effect of glycolysis, measure the blood sample as soon as possible after taking the sample. Samples of blood collected in vials with recommended anticoagulants should be analyzed within 30 minutes. Mix all anticoagulated samples thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8-10 times by hand before measurement. Alternatively follow local recommendations.

In cases of severe hypotension and peripheral circulatory failure, glucose measurement from capillary samples may be misleading. In such circumstances it is recommended to measure the glucose level using venous or arterial whole blood.⁶

Procedure

Read the Operating Manual¹ for the HemoCue Glucose 201 Analyzer for proper use of the analyzer.

Materials provided

HemoCue Glucose 201 Microcuvettes

Materials required but not provided

HemoCue Glucose 201 Analyzer

Lint-free tissue (non-fraying)

Lancet (capillary sampling)

Pipette or other transfer device (control material, venous or arterial samples)

Hydrophobic surface (control material, venous or arterial samples)

Commercial liquid quality controls recommended by HemoCue (contact HemoCue America Technical Support for advice)

Proper temperatures

Operating temperature for HemoCue Glucose 201 System and HemoCue Glucose 201 DM System is 64–86 °F (18–30 °C).

Directions for use

1. Pull the cuvette holder out to its loading position. Press and hold the on/off button until the display is activated. Read the relevant Operating Manual¹ for the analyzer for proper use of the system.
2. Use of the HemoCue Glucose 201 Microcuvettes.
 - a) The shape of the microcuvette and the names of the different parts of the microcuvette can be seen in figure 1.
 - b) Take the microcuvette out of the package. Hold the microcuvette opposite the filling end and bring into contact with the capillary specimen, see figure 2. When test tubes or syringes are used, the blood should be well mixed prior to the measurement. Place a drop of blood onto a hydrophobic surface with a transfer device, see figure 3. Always avoid contamination of the optical eye.
 - c) Fill the microcuvette completely by capillary action. **Do not refill the microcuvette!**
 - d) Wipe off the outside of the microcuvette with a clean lint free tissue, see figure 4. Do not touch the slit of the microcuvette. If air bubbles are present in the optical eye of a filled microcuvette, discard and fill a new microcuvette.
 - e) Place the microcuvette into the cuvette holder and start measurement as soon as possible but no later than 40 seconds after filling the microcuvette by gently pushing the cuvette holder to its measuring position, see figure 5.
 - f) After 40–240 seconds the HemoCue Glucose 201 Analyzer will find the steady-state of the chemical reaction and the result will appear on the display. The result will remain on the display as long as the cuvette holder is in the measuring position.
 - g) **Do not remeasure the microcuvette!**
 - h) Pull the cuvette holder out to its loading position and discard the used microcuvette. The microcuvette is for single use only.
 - i) If the display shows an error code an error has occurred, see the Troubleshooting Guide in the relevant Operating Manual¹ for the analyzer for more information.


Fig. 1


Fig. 2


Fig. 3

Fig. 4

Fig. 5

Quality control

The HemoCue Glucose 201 Analyzer has an internal electronic self-test. Every time the analyzer is turned on, it will automatically verify the performance of the optronic unit of the analyzer. This test is performed at regular intervals if the analyzer is left turned on.

The HemoCue Glucose 201 system must be verified on the days of testing using at least one level of commercially available controls. Follow local guidelines regarding quality control procedures. Only use controls recommended by HemoCue. For further information regarding controls contact HemoCue America/AB or your local distributor.

Results

The measured blood glucose value is read directly from the HemoCue Glucose 201 Analyzer in mg/dL and mmol/L. No calculations are necessary. The test is linear up to 444 mg/dL (24.6 mmol/L). Results above 444 mg/dL (24.6 mmol/L) are displayed as "HHH" or "overrange", and should be confirmed with another laboratory method.

To establish HemoCue glucose reference values and an intervention level, neonatal blood samples should be evaluated against a suitable laboratory method, taking into consideration the difference between whole blood and plasma reference values. Neonatal blood glucose levels below 50 mg/dL (2.8 mmol/L) in the first hours of life are cause for clinical concern and appropriate care should be exercised. Any results with HemoCue Glucose systems suggesting clinical intervention in the hyperglycemic range on pre-term neonates (<37 weeks), should be verified against a suitable laboratory reference method.

The following substances have been tested without interfering with the system: Ascorbic Acid*, Acetone*, Acetylsalicylic Acid**, Bilirubin*, Creatinine**, Urea**, Acetaminophen (Paracetamol)**, Dopamine**, Ephedrine**, Ibuprofen**, L-Dopa**, Methyl dopa**, Tetracycline**, Tolazamide**, Tolbutamide** and Uric Acid**.

* CLSI EP7-A2⁷

** CLSI EP7-P⁷

Limitations of the procedure

- The measurement needs to be started no later than 40 seconds after filling the microcuvette.
- Blood glucose results of the following specimens should be interpreted with caution: lipemic specimen (triglyceride, cholesterol and intralipid), toxic levels of methemoglobin (>5%) or high glucosamine concentration (>90 mg/dL, >5 mmol/L).
- Analyzers displaying plasma equivalent result multiply the measured whole blood glucose value by a factor of 1.11². Since this factor is based on normal hematocrit, care should be taken when the hematocrit may be extreme.

Expected values

Fasting glucose values (reference interval)⁸: Plasma glucose, adults 74-99 mg/dL (4.1-5.5 mmol/L).

For diagnosis of diabetes mellitus, follow local recommendations or use the following value according to WHO⁹:

Fasting plasma glucose, capillary or venous ≥ 126 mg/dL (≥ 7.0 mmol/L).

Specific performance characteristics

The results given below in "Within-run and total precision" come from one batch of microcuvettes and ten HemoCue Glucose 201 Analyzers. No restandardization was done during the analyzing period.

The results given in the "Correlation studies" come from HemoCue Analyzers which have not been restandardized during the study period. Studies A–D are whole blood glucose comparisons. Studies E and F are comparisons between the HemoCue Glucose 201 Analyzer with plasma conversion and a plasma glucose laboratory method. The neonatal studies (G–L) are whole blood glucose comparisons that have been performed using the HemoCue B-Glucose Analyzer. The HemoCue B-Glucose Analyzer is a system which is substantially equivalent to the HemoCue Glucose 201 system.

Within-run and total precision

Within-run precision and total precision were determined according to the CLSI standard EP5-A¹⁰. Commercially available freeze dried plasma controls at three different levels were used. The glucose concentration was measured in duplicate twice a day, morning and afternoon, during twenty consecutive days on all ten analyzers. Each day a new vial of each control level was used.

Plasma Control level	Number of determinations	Grand mean mg/dL	Grand mean mmol/dL	Within-run precision			Total precision		
				SD mg/dL	SD mmol/L	CV %	SD mg/dL	SD mmol/L	CV %
1	800	107.8	6.0	2.02	0.11	1.9	2.29	0.13	2.1
2	800	180.5	10.0	1.87	0.10	1.0	2.44	0.14	1.4
3	800	269.4	15.0	2.78	0.15	1.0	3.30	0.18	1.2

Correlations studies

The results of the comparison studies between HemoCue glucose and other glucose methods are summarized in the table below.

Adult studies	No of samples	Min. mg/dL	Min. mmol/L	Max. mg/dL	Max. mmol/L	Regression line mg/dL	Regression line mmol/L	Corr. coeff.
A	84	28	1.6	347	19.3	Y=1.04X+2	Y=1.04X+0.11	0.988
B	100	4	0.2	342	19.0	Y=1.01X+2	Y=1.01X+0.11	0.989
C	49	59	3.3	322	17.9	Y=1.06X-11	Y=1.06X-0.61	0.989
D	40	62	3.4	293	16.3	Y=1.00X-4	Y=1.00X-0.22	0.983
E	83	81	4.5	275	15.3	Y=1.05X-9	Y=1.05X-0.50	0.979
F	45	79	4.4	238	13.2	Y=0.99X-6	Y=0.99X-0.33	0.971

Neonatal studies	No of samples	Min. mg/dL	Min. mmol/L	Max. mg/dL	Max. mmol/L	Regression line mg/dL	Regression line mmol/L	Corr. coeff.
G	95	2	0.1	294	16.3	$Y=1.11X-6$	$Y=1.11X-0.33$	0.980
H	104	6	0.3	358	19.9	$Y=1.10X+2$	$Y=1.10X+0.11$	0.980
I	27	2	0.1	78	4.3	$Y=1.06X+3$	$Y=1.06X+0.17$	0.964
J	45	17	0.9	70	3.9	$Y=1.06X-5$	$Y=1.06X-0.28$	0.975
K	41	16	0.9	70	3.9	$Y=1.08X-5$	$Y=1.08X-0.28$	0.975
L	24	7	0.4	49	2.7	$Y=1.07X-4$	$Y=1.07X-0.22$	0.970

Note: Minimum (Min.) and maximum (Max.) values are taken from the comparative method. Also note that the neonatal samples contained hematocrit values between 30 % and 78 %.

References

- HemoCue Glucose 201 Operating Manual, HemoCue Glucose 201 DM Operating Manual.
- Fogh-Andersen et al, Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JIFCC, Vol 12 No 4.
- Burrin J.M., Alberti K.G.M.M, Diabetic Medicine, 1990; 7:199-206.
- Banauch et al, Z. Klin. Chem. u. Klin. Biochem, 1975; 13:101-107.
- Bergmeyer, Methods of Enzymatic Analysis 1, 1974; Chemie Publishers, Weinheim.
- Atkin et al, Annals of Internal Medicine, 1991, 114:12, 1020-1024.
- CLSI document EP7, Interference Testing in Clinical Chemistry.
- Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. 8th Edition; edited by Rifai, Horvath and Wittwer; 2018.
- Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycaemia. Report of a WHO/IDF Consultation. Geneva: WHO 2006.
- CLSI Document EP5-A, Vol 19, No 2, Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline.

Symbols used

Caution

Batch code

Use-by date
(Year Month Day)

In vitro diagnostic
medical device

Catalogue number

Rx only Prescription device
(applicable for US only)

Do not re-use

Temperature limitation

Date of manufacture

Consult instructions for
use

Once opened – Expiry date

Open vial expiry date. Must not
exceed the "Use-by" date

Manufacturer

Manufacturer

HemoCue AB
Kuvettgatan 1
SE-262 71 Ängelholm
Sweden

Phone +46 77 570 02 10
Fax +46 77 570 02 12

hemocue.com
info@hemocue.se

HemoCue Distributor USA

HemoCue America
250 South Kraemer Boulevard
Brea, CA, 92821
USA

Phone (general) 800.881.1611
Orders 800.323.1674
Technical support 800.426.7256
Fax (cust. service) 800.333.7043

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Références

- Manuel d'utilisation HemoCue Glucose 201, Manuel d'utilisation HemoCue Glucose 201 DM.
- Fogh-Andersen et al., Recommendations on reporting results for blood glucose (From an IFCC Stage 1 Document) IFCC Scientific Division Working Group on Selective Electrodes, JifCC, Vol 12 No 4.
- Burrin J.M., Alberti K.G.M.M, Diabetic Medicine, 1990; 7:199-206.
- Banauch et al., Z. Klin. Chem. u. Klin. Biochem, 1975; 13:101-107.
- Bergmeyer, Methods of Enzymatic Analysis 1, 1974; Chemie Publishers, Weinheim.
- Atkin et al., Annals of Internal Medicine, 1991, 114;12, 1020-1024.
- CLSI document EP7, Interference Testing in Clinical Chemistry.
- Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. 8th Edition; edited by Rifai, Horvath and Wittwer; 2018.
- Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycaemia. Report of a WHO/IDF Consultation. Geneva : WHO 2006.
- CLSI Document EP5-A, Vol 19, No 2, Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline.

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